



SPECIFICATION

PRODUCT	CODE	DIAMETER	LENGTH
	SITA03140	3mm	140mm
	SITA03330	3mm	330mm
	SITA05140	5mm	140mm
	SITA05280	5mm	280mm
	SITA05330	5mm	330mm
	SITA05380	5mm	380mm
	SITA05450	5mm	450mm
	SITA10140	10mm	140mm
	SITA10280	10mm	280mm
	SITA10330	10mm	330mm
	SITA10380	10mm	380mm
	SITA10450	10mm	450mm

DISPOSABLE SUCTION IRRIGATOR

Surgsci disposable suction irrigator features an excellent ergonomic design, allowing for multiple gripping methods and providing a more flexible operating experience.



Improves field visibility




Locates bleeding points




Reduces infection probability




DISPOSABLE SUCTION IRRIGATOR



STAINLESS STEEL TUBE




AVAILABLE IN VARIOUS SPECIFICATIONS

- A**  Product Diameter: $\phi 10$ mm
- B**  Product Diameter: $\phi 5$ mm
- C**  Product Diameter: $\phi 3$ mm

AVAILABLE IN VARIOUS SPECIFICATIONS:


Diameter: $\phi 3$ mm, $\phi 5$ mm, $\phi 10$ mm Length: 140mm, 280mm, 330mm, 380mm, 450mm

MULTIPLE HOLDING METHODS



Excellent ergonomic design provides surgeons with multiple flexible gripping options. The curved shape and anti-slip stripes optimize the experience of single-handed holding and operation. The sliding button is simple and easy to use, with a fixed position design that reduces fatigue during long-term operation, providing a good user experience.

USAGE SCENARIO



- Comprehensive Functionality**
Strong flushing and suction with blunt dissection and tissue manipulation, reducing instrument changes and enhancing surgical smoothness.
- Variety of Sizes**
Surgsci disposable suction irrigators are available in multiple sizes. The irrigation tube assembly supports replacement (future development includes components for electrocoagulation, electrotomy, and electric hooks).
- High Efficiency**
Independently designed internal sealing structure provides stable and reliable sealing. Fixed gear design reduces user fatigue.