


- Length 120mm /150mm
- Stainless steel
- With an internal spring-loaded blunt stylet
- having a coloured safety band
- Transparent handle
- with stopcock
- Luer lock connection
- Sterile
- Disposable



Visible Coloured safety band in the handle to indicate when the stylet mechanism is retracted for insertion and when it has automatically popped forward after peritoneal retraction

Supplied individually packed in a peel away pouch

| PRODUCT CATALOG | PRODUCT CODE | DESCRIPTION |
|-----------------|--------------|-------------|
| | VN120 | 120mm |
| | VN150 | 150mm |

Company address:

Area A, 4/F, Summer Valley, Mei Sheng Hui Gu Technology Park, No. 83 Dabao Road, Shanghe Community, Xin 'an Street, Bao 'an District, Shenzhen.

Factory Address :

3/F, Building 7, No.14 Zhongxing Road, Kengzi, Pingshan District, Shenzhen, Guangdong, China.

Tel: 0755-23707593

Email: service@surgsci.com